


Anterolateral Ligament Reconstruction

Surgical Technique



Anterolateral Ligament Reconstruction

In 1879, the French surgeon Paul Segond described a remarkably constant avulsion fracture pattern at the proximal-lateral tibia as a result of forced internal rotation. He described the existence of “a pearly, resistant, fibrous band” connecting the femur with the lateral tibia that seemed to show “extreme amounts of tension during forced internal rotation of the knee”. Later, Hughston described the importance of the “middle third capsular ligament” and how it was frequently torn in combination with an ACL tear, and this structure was essentially in this same anatomic location Segond described. However, only recently related to Claes’ 2013 publication, this anatomical structure has been comprehensively characterized as the anterolateral ligament (ALL).¹

Given its location at the anterolateral aspect of the knee, the ALL was found to act as an important restraint to internal tibial rotation. With the pivot-shift consisting of a coupled translation/rotation phenomenon, experimental sectioning of the ALL was found to invariably induce high-grade pivot-shifts in ACL-deficient cadaveric knees, unlike isolated ACL injury. In other words: high-grade pivot-shifts were only seen in the combination of ACL plus ALL injuries. With the continued occurrence of unsettled rotational laxity despite appropriately performed anatomic ACL reconstruction being a significant issue in current practice, the aforementioned new insights of the anatomy and function of the ALL could open the door to a potential solution. Specifically, reconstruction of the ALL could play a major role in improving results of isolated ACL reconstruction by providing better rotational control of the knee.

The goal of ALL reconstruction is to eliminate any residual rotational laxity and also reduce the risk of ACL graft rupture. Although contemporary ACL reconstruction is generally thought to deliver good results with excellent control of anterior-posterior (AP) laxity, the persistence of some degree of rotational instability characterized by a positive pivot-shift test in some patients is not uncommon. Rotational laxity after ACL injury is best quantified in the pivot-shift phenomenon, which is the most specific test that correlates best with functional outcome after reconstruction. However, the problematic persistence of a positive pivot-shift remains an unsolved issue in a significant amount of cases after both single and double ACL reconstruction.²

Historically, anterior laxity in anterior cruciate ligament deficient knees was treated surgically by isolated extraarticular iliotibial band tenodesis, as described by Lemair, Jacob and MacIntosh.^{3,4} This procedure in isolation was largely abandoned when arthroscopic single-bundle intraarticular ACL reconstruction emerged as the gold standard surgical treatment of ACL tear. More recently, the anterolateral ligament (ALL) has been shown to have an effect on rotational stability in several studies when performed in association with a standard intraarticular reconstruction of the ACL.^{5,6,7}

Anterolateral Ligament Reconstruction Indications

Anterolateral ligament reconstruction is aimed at augmenting rotational stability in the ACL reconstructed knee. Because combined injuries to both the ACL, and ALL or deep IT Band, act as a prerequisite for the occurrence of an IKDC grade III pivot-shift, ACL-injured patients with a high-grade pivot-shift might benefit from an additional anterolateral reconstruction in order to avoid persistent rotational laxity. Hyperlax females with excessive recurvatum and physiologic joint laxity are potentially appropriate candidates for combined ACL reconstruction and extra-articular stabilization. Furthermore, in ACL-injured pivoting athletes who require absolute stability, anterolateral reconstruction should be contemplated if only an IKDC grade II pivot-shift is present. Finally, revision ALR reconstruction cases commonly exhibit significant rotational laxity due to a tendency for increased joint laxity from previous meniscus removal or resultant laxity of secondary ligamentous restraints. Especially in the absence of frank re-trauma or obvious technical errors explaining graft failure, concomitant ALR reconstruction should always be considered as a means of improving stability in these complex cases.

¹ Claes et al., *Journal of Anatomy*, 2013.

² Suomalainen et al., *AJSM*, 2012.

³ M. Lemaire et al., *Rev Chir Orthop Reparatrice Appar Mot*, 1980.

⁴ J. Ireland et al., *J Bone Joint Surg Br*, 1980.


⁵ Sonnery-Cottet et al., *American Journal of Sports Medicine*, 2015.

⁶ Pomajzl et al., *Journal of Arthroscopic and Related Surgery*, 2014.

⁷ Parsons et al., *American Journal of Sports Medicine*, 2015.

⁸ Helito et al., *The American Journal of Sports Medicine*, 2014.

⁹ Amis et al., *The Journal of Bone and Joint Surgery*, 2013.



FCL femoral attachment

ALL femoral attachment

Popliteus tendon

LE

GERDY TUBERCLE

Biceps femoris


ALL tibial attachment

Lateral gastrocnemius

Anterior arm of biceps femoris


FCL attachment to fibular head

Anterolateral Ligament Reconstruction




The anterolateral ligament (ALL) will run obliquely and in a double stranded construct to the anterolateral tibia, where it is transosseous wise looped and lead back to the femoral insertion posterior and proximal to the epicondylus. The femoral stab incision is made slightly proximal and posterior to the lateral epicondyle. The two tibial stab incisions are made proximal to the anterior edge of the fibular head and 18 mm anterior to this incision.


A gracilis tendon autograft is utilized, as the cross-sectional area, length, and ultimate loads provide sufficient anatomical and mechanical properties for ALL reconstructions.⁷ Minimum graft length is 24 cm, whipstitched 20 mm at the thin, proximal end with a 2-0 FiberWire suture. The graft should be tapered at the whipstitched end to facilitate insertion of the graft into the femur.



On the femur, the insertion point for the 2.4 mm guide pin is about 8 mm proximal and 4 mm posterior to the lateral epicondyle. When drilling the 2.4 mm pin, aim slightly anterior and proximal in order to avoid drilling into the femoral socket of the ACL reconstruction.^{8,9}




After the pin is drilled, split the iliotibial band around it with a scalpel to facilitate the insertion of a SwiveLock anchor.




4

Use the 18 mm parallel drill guide and place the first 2.4 mm guide pin through the posterior stab incision, 10 mm distal to the joint line at the anterior aspect of the fibular head . Place the second 2.4 mm drill pin through the anterior sleeve, being parallel to the tibial joint line. The anterior pin should be placed between the center and the posterior margin of the Gerdy's tubercle.




5


The tensioning behaviour of the graft can be tested by wrapping a FiberWire around the pins in full extension and moving the leg through the range-of-motion. When the FiberWire is held taut at extension, it should be isometric or become slightly lax by no more than 5 mm as it moves to 90° of flexion.





6

Overdrill the guide pins in the tibia with a 5 mm cannulated drill to a depth of 15 mm.






The 90° transosseous obturator (AR-5524) is used to connect the anterior with the posterior tibial socket. To protect the surrounding bone, the max leverage of the obturator is 90° to the bone. The same is done from posterior to anterior.




The transosseous obturator is used to pass a shuttle suture from the posterior to the anterior tibial hole.

The suture is captured with a forceps.



10


Overdrill the femoral guide pin with a 5 mm cannulated drill to a depth of 20 mm. For hard bone, the 4.75 mm SwiveLock tap can be used. Place a nitinol guide wire in the femoral socket.



11


Load the whipstitched sutures of the gracilis graft through the PEEK eyelet of the 4.75 mm SwiveLock anchor.

Insert the SwiveLock with the gracilis graft end into the drilled hole. Ensure the eyelet is fully seated so that the threads of the SwiveLock anchor have started to touch the bone. Hold the paddle of the SwiveLock and turn the knob to insert the anchor. The SwiveLock sutures can be removed.




12

Use a curved hemostat to dissect underneath the iliotibial band to create a plane from the femoral incision to the tibial incisions. Place a passing suture from distal to proximal with the hemostat. Use the passing suture to pass the graft to the tibial side.




13


Use the tibial transosseous shuttle suture to pass the graft through the osseous tunnel from the posterior to the anterior hole. For screw insertion place a nitinol wire in the anterior tibial socket.




Using a femoral passing suture, the graft is passed again proximal, underneath the iliotibial band towards the femoral insertion.



With the knee in full extension and neutral rotation, the graft is fixed in the anterior tibial socket using a 4.75 mm × 15 mm BioComposite tenodesis screw.



Again with the knee in full extension and neutral rotation, the proximal end of the graft is fixed to the femoral insertion using the FiberWire suture coming from the SwiveLock cannulation.



Ordering Information

Anterolateral Ligament Reconstruction

Product Description	Item Number
Transosseous drill guide, 18 mm	AR-5523
Transosseous obturator for Ø 5 mm	AR-5524
Drill tip guide pin, 2.4 mm × 311 mm	AR-1250L-1
Cannulated drill, 5 mm	AR-1205L
Nitinol guide pin for bio-interference screws, 1.1 mm	AR-1249
BioComposite SwiveLock C vented, with / closed eyelet, 4.75 mm × 19.1 mm	AR-2324BCC
BioComposite tenodesis screw, 4.75 mm × 15 mm	AR-1547BC
FiberLoop® # 0, with straight needle	AR-7253
# 2 FiberWire suture	AR-7233

Products advertised in this brochure/surgical technique guide may not be available in all countries. For information on availability, please contact Arthrex Customer Service or your local Arthrex representative.

*Data on file



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's Directions For Use. Postoperative management is patient specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level and/or outcomes.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

© Arthrex GmbH, 2019. All rights reserved. | www.arthrex.com | LT2-80003-en-US_A